

Watt lab

Single cell cDNA library protocol

Prepare beforehand

1. 1mM dNTP mix from 100mM dATP, dGTP, dCTP, dTTP (separately from Thermo Scientific)
2. Spike RNA mix from *A. thaliana* mRNA LTP4, LTP6, NAC1, TIM (separately from Stratagene) diluted to 10^{-1} pg/ μ l, 10^{-2} pg/ μ l, 10^{-3} pg/ μ l, 10^{-4} pg/ μ l, respectively

First-strand synthesis

First-strand synthesis buffer:

5x SuperScript III buffer (Invitrogen)	1 μ l
5% NP-40 (Pierce)	0.5 μ l
1mM dNTP mix (Thermo Scientific)	0.25 μ l
1 μ M MO ₄ d(T) primer*	0.075 μ l
0.1M DTT (Invitrogen)	0.05 μ l
SuperRNaseIN (Ambion)	0.25 μ l
RNaseOUT (Invitrogen)	0.25 μ l
Spike RNA mix (Stratagene)	1.0 μ l
Nuclease-free water (Ambion)	0.625 μ l
Total	<hr/> 4 μ l

*MO₄d(T) primer sequence (Osawa et al., Development, 2005):

AAG CAG TGG TAT CAA CGC AGA GTG GCC ATT ACG GCC GTA CTT TTT
TTT TTT TTT TTT TTT TTT TTT TTT T (DNA Technology, HPLC-purified)

1. Aliquot 4 μ l of first-strand synthesis buffer into thin-wall PCR tubes on ice.
2. Seed single cell in 0.5 μ l into tube.
3. Freeze immediately on dry ice.
4. Incubate at 65°C for 5 min to lyse cell.
5. Incubate at 45°C for 2 min to allow the primer to anneal.
6. Add 0.5 μ l SuperScript III (Invitrogen).
7. Incubate at 45°C for 15 min.
8. Incubate at 65°C for 10 min to inactivate the RT enzyme.

Exonuclease I treatment

Exonuclease I reaction mix:

Exonuclease I (Thermo Scientific)	0.2 μ l
75mM MgCl ₂ (Ambion)	0.6 μ l
Nuclease-free water (Ambion)	0.2 μ l
Total	<hr/> 1 μ l

1. Add 1 μ l to each tube.
2. Incubate at 37°C for 30 min.
3. Incubate at 80°C for 25 min to inactivate the enzyme.

RNase H treatment and poly-adenylation

Reaction mix:

5x TdT buffer (Promega)	2.6 μ l
100mM dATP (Thermo Scientific)	0.097 μ l
RNase H (Invitrogen)	0.5 μ l
TdT (Promega)	0.5 μ l
Nuclease-free water (Ambion)	3.303 μ l
Total	<hr/> 7 μ l

1. Add 7 μ l to each tube.
2. Incubate at 37°C for 15 min.
3. Incubate at 70°C for 10 min to inactivate the enzymes.

PCR amplification 1

PCR reaction mix:

10x ExTaq buffer (Takara)	2 μ l
100mM dATP (Thermo Scientific)	0.13 μ l
100mM dGTP (Thermo Scientific)	0.13 μ l
100mM dCTP (Thermo Scientific)	0.13 μ l
100mM dTTP (Thermo Scientific)	0.13 μ l
100 μ M MO ₄ d(T) primer	1.65 μ l
ExTaq polymerase (Takara)	0.2 μ l
Nuclease-free water (Ambion)	11.63 μ l
Total	<hr/> 16 μ l

1. Aliquot 16 μ l into fresh tubes; add 4 μ l poly-adenylated cDNA in triplicate.

2. Perform PCR:

94°C	1 min		1 cycle
50°C	2 min		
72°C	2 min		
94°C	30 sec		35 cycles
60°C	30 sec		
72°C	2 min		

3. Pool the triplicate PCR products.

PCR amplification 2

PCR reaction mix:

10x ExTaq buffer (Takara)	2µl
2.5mM dNTP mix (Takara)	1.6µl
100µM T7-MO ₄ primer*	0.4µl
ExTaq polymerase (Takara)	0.2µl
Nuclease-free water (Ambion)	13.8µl
Total	<hr/> 18µl

*T7-MO4 primer:

GGCCAGTGAATTGTAATACGACTCACTATAGGGAGGCGGAAGCAGTGGTATC
AACGCAGAGTGGCCATTACGGCCGTAC (DNA Technology, HPLC-purified)

1. Aliquot 18µl into fresh tubes; add 2µl pooled PCR 1 product in duplicate.
2. Perform PCR:

94°C	30 sec		35 cycles
60°C	30 sec		
72°C	2 min		

3. Pool the duplicate PCR products.

PCR product purification

Purify PCR product with QIAquick PCR Purification Kit (QIAGEN).

Analyse amplified cDNA

Run 5µl PCR product + 1µl 6x DNA loading buffer on a 2% agarose gel. A smear 500-700bp in size should be seen. The amplified cDNA can now be used for the in vitro transcription reaction to generate labelled cRNA for hybridisation to microarrays.